

Product Datasheet

Anti-IL-3Ra / CD123 Reference Antibody (orb1806264)

Description

Anti-IL-3Ra / CD123 Reference Antibody is

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

IL-3Ra / CD123

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

146.06 kDa


Expiration Date

12 months from date of receipt.


IL-3Ra NR M R



Anti-IL-3Ra
/ CD123
Reference
Antibody
(...



The purity
of Anti-IL-
3Ra /
CD123
Refere...



Immobilized
human IL
3Rα His at
2 µg/m...