



Product Datasheet

Anti-IL-3Ra / CD123 Reference Antibody (orb1806265)

Description	Anti-IL-3Ra / CD123 Reference Antibody is expressed from CHO. It has a predicted MW of 146.04 kDa.
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	IL-3Ra / CD123
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	146.04 kDa
Expiration Date	12 months from date of receipt.



Anti-IL-3Ra
/ CD123
Reference
Antibody
(...



The purity
of Anti-IL-
3Ra /
CD123
Refere...