

Product Datasheet

Anti-CCR2 / CD192 Reference Antibody (orb1806269)

Description Anti-CCR2 / CD192 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the ...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CCR2 / CD192

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 150 kDa

Expiration Date 12 months from date of receipt.



Anti-CCR2 /
CD192
Reference
Antibody
(pl...



The purity
of Anti-
CCR2 /
CD192
Referenc...