



Product Datasheet

Anti-C1s Reference Antibody (orb1806275)


Description	Anti-C1s Reference Antibody is expressed from CHO. The heavy chain type is hulg4SP, and the light c...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	C1s
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.



Anti-C1s
Reference
Antibody
(sutimlimab)...



The purity of
Anti-C1s
Reference
Antibod...



Immobilized
human C1s
His at 2 µg/mL
ca...