

## Product Datasheet

### Anti-B7-H6 / NCR3LG1 Reference Antibody (orb1806281)

**Description** Anti-B7-H6 / NCR3LG1 Reference Antibody is expressed from CHO. It has a predicted MW of 145.5 kDa.

**Reactivity** Human

**Endotoxins** < 1 EU/mg

**Conjugation** Unconjugated

**Tested Applications** ELISA, FA, FACS, Kinetics

**Target** B7-H6 / NCR3LG1

**Preservatives** 100 mM Pro-Ac 20 mM Arg pH 5.0

**Storage** -80°C

**Note** For research use only


**Isotype** IgG

**Clonality** Monoclonal


**Purity** > 95%

**MW** 150 kDa

**Expiration Date** 12 months from date of receipt.



Anti-B7-H6  
/ NCR3LG1  
Reference  
Antibody ...



The purity  
of Anti-B7-  
H6 /  
NCR3LG1  
Refer...