


Product Datasheet

Human Tetranectin/CLEC3B ELISA Kit (orb180629)

Description	Human Tetranectin/CLEC3B ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E22-V202
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants and serum.
MW	22537 MW
Uniprot ID	P05452
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:2000, concentration in serum is 2-10 ug/ml.
Expiration Date	Please enquire.



Human
Tetranectin/CLEC3B
PicoKine ELISA ...