

Product Datasheet


Anti-TNFRSF9 / 4-1BB / CD137 Reference Antibody (orb1806297)

Description

Anti-TNFRSF9 / 4-1BB / CD137 Reference

	chain type is hulg...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFRSF9 / 4-1BB / CD137
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	146.5 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-TNFRSF9 / 4-1BB / CD137 Reference A...



The purity of Anti-TNFRSF9 / 4-1BB / CD1...



Immobilized N His HSA TNFSF9/4 1BB at 2 ...