


Product Datasheet

Human Thioredoxin ELISA Kit (orb180630)

Description	Human Thioredoxin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: V2-V105
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11737 MW
Uniprot ID	P10599
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human
Thioredoxin
PicoKine
ELISA Kit
sta...