


Product Datasheet

Human CD320 / 8D6A / TCbIR ELISA Kit (orb180631)

Description	Human CD320 / 8D6A / TCbIR ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S36-Y229
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	28991 MW
Uniprot ID	Q9NPF0
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 100-360 pg/ml.
Expiration Date	Please enquire.



Human
CD320/8D6A
PicoKine
ELISA Kit
stan...