


Product Datasheet

Human Ficolin-1 / FCN1 / M Ficolin ELISA Kit (orb180632)

Description	Human Ficolin-1 / FCN1 / M Ficolin ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A30-A326
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	35078 MW
Uniprot ID	000602
Sensitivity	20 pg/ml
Dilution Range	Dilution ratio of 1:20, concentration in serum and plasma is 120 ng/ml.
Expiration Date	Please enquire.



Human
Ficolin-1
PicoKine
ELISA Kit
stand...