


Product Datasheet

Human Ficolin-2 ELISA Kit (orb180633)

Description	Human Ficolin-2 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: L26-A313
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	34001 MW
Uniprot ID	Q15485
Sensitivity	20 pg/ml
Expiration Date	Please enquire.



Human
Ficolin-2
PicoKine
ELISA Kit
stand...