


Product Datasheet

Human Ficolin-3 ELISA Kit (orb180634)

Description	Human Ficolin-3 ELISA Kit
Reactivity	Human
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q24-R299
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	32903 MW
Uniprot ID	075636
Sensitivity	20 pg/ml
Dilution Range	For serum, dilution ratio is 1:1 and concentration is around 7 ng/ml. For plasma, dilution ratio is 1:5 and concentration is around 15 ng/ml.
Expiration Date	Please enquire.



Human
Ficolin-3
PicoKine
ELISA Kit
stand...