


Product Datasheet

Anti-ERBB2 / HER2 / CD340 Reference Antibody (orb1806343)

Description

Anti-ERBB2 / HER2 / CD340 Reference Antibody is


	145.5 k...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	ERBB2 / HER2 / CD340
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.



Anti-ERBB2
/ HER2 /
CD340
Reference
Anti...



The purity
of Anti-
ERBB2 /
HER2 /
CD340 ...



Immobilized
human
HER2; His
at 2 µg/mL
...