


Product Datasheet

Anti-GPC3 / Glypican-3 Reference Antibody (orb1806346)

Description

Anti-GPC3 / Glypican-3 Reference Antibody is


	hulgG1, and...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	GPC3 / Glypican-3
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.



Anti-GPC3 /
Glypican-3
Reference
Antibod...



The purity
of Anti-
GPC3 /
Glypican-3
Ref...



Immobilized
human
GPC3; His
tag at 2
µg...