



Product Datasheet

Anti-FOLR1 Reference Antibody (orb1806347)


Description	Anti-FOLR1 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light C...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	FOLR1
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.2 kDa
Expiration Date	12 months from date of receipt.



Anti-FOLR1 /
FRA
Reference
Antibody
(mir...



The purity
of Anti-
FOLR1 / FRA
Reference...



Immobilized
human FRa
His at 2
µg/mL ca...