

Product Datasheet


Anti-FGFR2 / CD332 Reference Antibody (orb1806349)

Description

Anti-FGFR2 / CD332 Reference Antibody is

	hulG1, and the...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	FGFR2 / CD332
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	78.2 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-FGFR2
/ CD332
Reference
Antibody
(a...



The purity
of Anti-
FGFR2 /
CD332
Referen...



Human
FGFR2-β-IIIc
(A13) CHOS
cells wer...