


Product Datasheet

Human IL-18R1 ELISA Kit (orb180635)

Description	Human IL-18R1 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: C22-R329
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	62304 MW
Uniprot ID	Q13478
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-18R1 PicoKine ELISA Kit standard...