

Product Datasheet


Anti-FGFR2 / CD332 Reference Antibody (orb1806350)

Description

Anti-FGFR2 / CD332 Reference Antibody is

	hulgG1, and the...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	FGFR2 / CD332
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	145.24 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-FGFR2 / CD332 Reference Antibody (b...



The purity of Anti-FGFR2 / CD332 Referen...



SNU-16 cells were stained with Anti-FGFR...