

Product Datasheet

Anti-CXCR4 / CD184 Reference Antibody (orb1806361)

Description

Anti-CXCR4 / CD184 Reference Antibody is

Reactivity

Human

Endotoxins

< 1 EU/mg

Conjugation

Unconjugated

Tested Applications

ELISA, FA, FACS, Kinetics

Target

CXCR4 / CD184

Preservatives

100 mM Pro-Ac 20 mM Arg pH 5.0

Storage

-80°C

Note

For research use only

Isotype

IgG

Clonality

Monoclonal

Purity

> 95%


MW

144.84 kDa

Expiration Date

12 months from date of receipt.


kDa NR M R



Anti-CXCR4
/ CD184
Reference
Antibody
(u...



The purity
of Anti-
CXCR4 /
CD184
Referen...



Human
CXCR4
Romas cells
were
stained
wit...