

Product Datasheet


Anti-TNFSF7 / CD27L / CD70 Reference Antibody (orb1806368)

Description


Anti-TNFSF7 / CD27L / CD70 Reference Antibody

	hulG1,...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFSF7 / CD27L / CD70
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.


kDa NR M R



Anti-TNFSF7 / CD27L / CD70 Reference Antibody



The purity of Anti-TNFSF7 / CD27L / CD70...



HL-60 cells were stained with Anti-TNFSF...