


Product Datasheet

Anti-TNFRSF7 / CD27 Reference Antibody (orb1806375)


Description

Anti-TNFRSF7 / CD27 Reference Antibody is


	hulgG1, and th...
Reactivity	Human
Endotoxins	< 1 EU/mg
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FACS, Kinetics
Target	TNFRSF7 / CD27
Preservatives	100 mM Pro-Ac 20 mM Arg pH 5.0
Storage	-80°C
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Purity	> 95%
MW	150 kDa
Expiration Date	12 months from date of receipt.



Anti-TNFRSF7 / CD27 Reference Antibody (...)



The purity of Anti-TNFRSF7 / CD27 Refere...



Immobilized human CD27 FC at 2 µg/mL ca...