

Product Datasheet

Anti-CD20 Reference Antibody (orb1806379)

Description Anti-CD20 Reference Antibody is expressed from CHO. The heavy chain type is hulgG1, and the light ch...

Reactivity Human

Endotoxins < 1 EU/mg

Conjugation Unconjugated

Tested Applications ELISA, FA, FACS, Kinetics

Target CD20

Preservatives 100 mM Pro-Ac 20 mM Arg pH 5.0

Storage -80°C

Note For research use only


Isotype IgG

Clonality Monoclonal


Purity > 95%

MW 146.32 kDa


Expiration Date 12 months from date of receipt.



Anti-CD20
Reference
Antibody
(rituximab)...



The purity of
Anti-CD20
Reference
Antibo...



Immobilized
human EGFR
His at 2
µg/mL c...