


Product Datasheet

Human Prostatin ELISA Kit (orb180638)

Description	Human Prostatin ELISA Kit
Reactivity	Human
Range	1.56 ng/ml - 100 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: A33-G319
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	36431 MW
Uniprot ID	Q16651
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 1 ng/ml.
Expiration Date	Please enquire.



Human
Prostatin
PicoKine
ELISA Kit
stand...