


Product Datasheet

Mouse Prostatin ELISA Kit (orb180639)

Description	Mouse Prostatin ELISA Kit
Reactivity	Mouse
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: A30-Q289
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	36729 MW
Uniprot ID	Q9ESD1
Sensitivity	50 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is less than the lowest standard, 0.78 ng/ml.
Expiration Date	Please enquire.



Mouse
Prostatin
PicoKine
ELISA Kit
stand...