



Product Datasheet

CD49f Antibody (APC) (orb1806413)

Description	CD49f (alpha 6 integrin) is a type I transmembrane glycoprotein, which noncovalently associates with...
Reactivity	Bovine, Canine, Equine, Feline, Human, Mouse, Porcine, Primate, Rabbit, Sheep
Conjugation	APC
Tested Applications	FC
Immunogen	mouse mammary tumor cells
Target	CD49f
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.
Isotype	Rat IgG2a
Clonality	Monoclonal
Uniprot ID	P23229
Entrez	3655
Expiration Date	12 months from date of receipt.



Separation of human CD49f positive throm...



Flow cytometry surface staining pattern ...