


Product Datasheet

Human DLL1 ELISA Kit (orb180642)

Description	Human DLL1 ELISA Kit
Reactivity	Human
Range	78 pg/ml - 5,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S22-G540
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , urine and human milk.
MW	78056 MW
Uniprot ID	000548
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum is 830-2162 pg/ml.
Expiration Date	Please enquire.



Human DLL1
PicoKine
ELISA Kit
standard c...