



## Product Datasheet

### IgM Antibody (orb1806429)


|                            |                                                                                                                                                                                                                                                     |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | IgM Antibody                                                                                                                                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                        |
| <b>Tested Applications</b> | FC, ICC, IHC, IP, WB                                                                                                                                                                                                                                |
| <b>Immunogen</b>           | Human IgM                                                                                                                                                                                                                                           |
| <b>Target</b>              | IgM                                                                                                                                                                                                                                                 |
| <b>Preservatives</b>       | Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide                                                                                                                                                                       |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                                                                                                           |
| <b>Concentration</b>       | 100 µg/ml                                                                                                                                                                                                                                           |
| <b>Storage</b>             | Shelf life: 1 year from date of receipt. Storage: 2 - 8°C                                                                                                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                               |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                                                                                                   |
| <b>Clonality</b>           | Recombinant                                                                                                                                                                                                                                         |
| <b>Dilution Range</b>      | WB: 1:1000, IP: 20 µl/mg lysate, IHC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections, ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                     |



Detection of human IgM by western blot. ...



Detection of mouse IgM by western blot. ...



Detection of human IgM by immunohistoche...