


Product Datasheet

Mouse DLL4 ELISA Kit (orb180643)

Description	Mouse DLL4 ELISA Kit
Reactivity	Mouse
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S28-P525
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	74991 MW
Uniprot ID	Q9J171
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 60-80 pg/ml.
Expiration Date	Please enquire.



Mouse DLL4
PicoKine
ELISA Kit
standard c...