



Product Datasheet

Cytokeratin 14 Antibody (orb1806432)


Description	Cytokeratin 14 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IP, WB
Immunogen	Human Cytokeratin 14 (C terminal)
Target	Cytokeratin 14
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P02533
Dilution Range	WB: 1:1000, IP: 20 µl/mg lysate, IHC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections, ICC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE cell sections
Expiration Date	12 months from date of receipt.



Detection of
human
Cytokeratin
14 by wes...



Detection of
human
Cytokeratin
14 by imm...



Detection of
human
Cytokeratin
14 (shade...