


Product Datasheet

Human Galectin-2 ELISA Kit (orb180644)

Description	Human Galectin-2 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: M1-E132
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	14644 MW
Uniprot ID	P05162
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 740-1226 pg/ml.
Expiration Date	Please enquire.



Human
Galectin-2
PicoKine
ELISA Kit
stan...