



Product Datasheet

CD48 Antibody (orb1806445)


Description	CD48 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	residues 29-220 (ECD)
Target	CD48
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P09326
Dilution Range	WB: 1:1000, IP: 20 µl/mg lysate
Expiration Date	12 months from date of receipt.



Detection of human CD48 by western blot....



Detection of human CD48 (shaded) in NK-92...



Detection of human CD48 by western blot ...