


Product Datasheet

Mouse Galectin-7/LGALS7 ELISA Kit (orb180645)

Description	Mouse Galectin-7/LGALS7 ELISA Kit
Reactivity	Mouse
Range	93.8 pg/ml - 6,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S2-F136
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	15173 MW
Uniprot ID	O54974
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
Galectin-7
PicoKine
ELISA Kit
stan...