



Product Datasheet

Troponin I Antibody (orb1806455)


Description	Troponin I Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	between 1 and 50
Target	Troponin I
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P19429
Dilution Range	WB: 1:1000, IHC: 1:100 to 1:500. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections
Expiration Date	12 months from date of receipt.



Detection of human Troponin I by western...



Detection of mouse Troponin I by western...



Detection of human Troponin I by immunoh...