


Product Datasheet

Human Corin ELISA Kit (orb180649)

Description	Human Corin ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: S67-N1042
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
MW	116486 MW
Uniprot ID	Q9Y5Q5
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 400-1000 pg/ml.
Expiration Date	Please enquire.



Human Corin
PicoKine
ELISA Kit
standard ...