

Product Datasheet

Mouse IL-6 R alpha / CD126 Protein, Fc Tag (MALS verified) (orb1806493)

Description

Mouse IL-6 R alpha Protein, Fc Tag is

Reactivity

Mouse

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

IL-6 R alpha

Preservatives

Powder

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 63.9 kDa. The protein migrates as 70-90 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

90%

Source

Mouse IL-6 R alpha Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 20 - Glu 357 (Accession # P22272-1).

MW


63.9 kDa

Uniprot ID

P22272

Expiration Date

6 months from date of receipt.



Immobilized
Mouse IL-6,
His Tag
(Cat. No...