


Product Datasheet

Human REG-4 ELISA Kit (orb180655)

Description	Human REG-4 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: D23-P158
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	18230 MW
Uniprot ID	Q9BYZ8
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 2-4 ng/ml.
Expiration Date	Please enquire.



Human REG-4 ELISA Kit standard curve