


Product Datasheet

Rat ADAM12 ELISA Kit (orb180657)

Description	Rat ADAM12 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: CHO; Immunogen sequence: R29-S513
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin).
MW	96518 MW
Uniprot ID	F1M9Z9
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is 320-520 pg/ml.
Expiration Date	Please enquire.



Rat ADAM12
PicoKine
ELISA Kit
standard c...