


Product Datasheet

Rat HGF/Hepatocyte growth factor ELISA Kit (orb180665)

Description	Rat HGF/Hepatocyte growth factor ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q33-R495 (alpha) & V496-L728 (beta)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA, citrate).
MW	82906 MW
Uniprot ID	P17945
Sensitivity	10 pg/ml
Dilution Range	For serum, dilution ratio is 1:20 and concentration is 20-40 ng/ml. For plasma, dilution ratio is 1:5 and concentration is 13192-50523 pg/ml.
Expiration Date	Please enquire.



Rat HGF
PicoKine
ELISA Kit
standard
curv...