


Product Datasheet

Human Cathepsin E ELISA Kit (orb180670)

Description	Human Cathepsin E ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q18-P396
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	43312 MW
Uniprot ID	P14091
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:1, concentration in serum and plasma is around 500 pg/ml.
Expiration Date	Please enquire.



Human
Cathepsin E
PicoKine
ELISA Kit
sta...