


Product Datasheet

Mouse TNFSF18/GITRL ELISA Kit (orb180672)

Description	Mouse TNFSF18/GITRL ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: T47-S173
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	19732 MW
Uniprot ID	Q7TS55
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Mouse
TNFSF18/GITRL
PicoKine ELISA
Kit s...