

Product Datasheet

Human Fc gamma RIIB/C (CD32b/c)(FCGR2B/C) ELISA Kit (orb180674)

Description Human Fc gamma RIIB/C (CD32b/c)(FCGR2B/C)
ELISA Kit

Reactivity Human

Range 46.9 pg/ml - 3,000 pg/ml

Tested Applications ELISA

Immunogen Expression system for standard: NS0;
Immunogen sequence: A46-P217

Storage Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)

Note For research use only

Sample Types cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).


MW 35578 MW

Uniprot ID **P31994, P31995**

Sensitivity 10 pg/ml

Dilution Range Dilution ratio of 1:1, concentration in serum and plasma is 562-1200 pg/ml.

Expiration Date Please enquire.



Human
CD32/FCGR2b/c
PicoKine ELISA
Kit s...