

## Product Datasheet

### Human CPP(Copeptin) ELISA Kit (orb1807117)

<b>Description</b>	Human CPP(Copeptin) ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	31.25-2000 pg/mL
<b>Target</b>	CPP
<b>Note</b>	For research use only
<b>Application notes</b>	<p>This ELISA kit uses the Competitive-ELISA principle. The micro ELISA plate provided in this kit has been pre-coated with Human CPP. During the reaction, Human CPP in the sample or standard competes with a fixed amount of Human CPP on the solid phase supporter for sites on the Biotinylated Detection Ab specific to Human CPP. Excess conjugate and unbound sample or standard are washed away, and Avidin-Horseradish Peroxidase (HRP) conjugate are added to each micro plate well and incubated. Then a TMB substrate solution is added to each well. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns from blue to yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm <math>\pm</math> 2 nm. The concentration of Human CPP in tested samples can be calculated by comparing the OD of the samples to the standard curve.</p>
<b>Sample Types</b>	serum, plasma, Tissue homogenate and Other biological samples
<b>Sensitivity</b>	13.7 pg/mL
<b>Expiration Date</b>	Please enquire.

