

Product Datasheet

Metalloproteinase inhibitor 2 TIMP2 Antibody (orb18073)

Description

Metalloproteinase inhibitor 2 TIMP2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Immunogen

A synthetic peptide corresponding to a sequence at the C-terminus of human TIMP2 (205-220aa RGAAPPKQEFLDIEDP), identical to the related mouse and rat sequences.

Western
Blot
analysis of
rat kidney
tiss...

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

24399 MW

Uniprot ID
P16035
Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Mouse, Rat, By Heat ELISA , 0.1-0.5µg/ml, Human

Expiration Date

12 months from date of receipt.