

Product Datasheet

Recombinant human FIBIN protein, N-Trx-His (orb1808083)

Description	Recombinant human FIBIN protein, N-Trx-His
Target	FIBIN
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	N-Trx-His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Protein Sequence	19-211/211
Purity	> 90% as determined by SDS-PAGE
MW	40.1 kDa
Uniprot ID	Q8TAL6
Expiration Date	6 months from date of receipt.