

## Product Datasheet

### General 8-OHdG ELISA Kit (orb1808938)

|                        |                                           |
|------------------------|-------------------------------------------|
| <b>Description</b>     | General 8-OHdG ELISA Kit                  |
| <b>Reactivity</b>      | All                                       |
| <b>Range</b>           | 62.5-4000pg/mL                            |
| <b>Storage</b>         | Store at 4°C & -20°C for 12 months.       |
| <b>Note</b>            | For research use only                     |
| <b>Sample Types</b>    | Serum, plasma or other biological fluids. |
| <b>Sensitivity</b>     | 18.5pg/mL                                 |
| <b>Expiration Date</b> | Please enquire.                           |