

Product Datasheet

TJP1 Antibody (orb1809436)

Description	TJP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ZO-1 (1361-1460/1748aa)
Target	TJP1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	191 kDa
Uniprot ID	Q07157
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.