


Product Datasheet

BI-78D3 (orb180962)

Description	BI-78D3 is a competitive c-Jun N-terminal kinase (JNK) inhibitor (IC ₅₀ = 280 nM) that displays > 100...
Note	For research use only
Formula	C ₁₃ H ₉ N ₅ O ₅ S ₂
MW	379.3711
CAS Number	[883065-90-5]
SMILES	<chem>O=C1NN=C(SC2=NC=C([N+])([O-])=O)S2)N1C3=CC=C(OCCO4)C4=C3</chem>
Expiration Date	12 months from date of receipt.



Chemical structure of BI-78D3