

Product Datasheet

Proteasome 20S beta 6 antibody (orb1810)

Description	Proteasome 20S beta 6 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Proteasome 20S beta 6 (21-120/239 aa)
Target	Proteasome 20S beta 6
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	P28072
Expiration Date	12 months from date of receipt.