

Product Datasheet

ICILIN (orb181068)

Description	Icilin is a synthetic super-agonist of the transient receptor potential M8 (TRPM8) ion channel.
Note	For research use only
Formula	$C_{16}H_{13}N_3O_4$
MW	311.2921
CAS Number	[36945-98-9]
SMILES	<chem>O=C1N([H])C(C2C([H])=C([H])C([H])=C(C=2[H])[N+](=O)[O-])=C([H])C([H])([H])N1C1=C([H])C([H])=C([H])C([H])=C1O[H]</chem>
Expiration Date	12 months from date of receipt.