

Product Datasheet

OSI 930 (orb181074)

Description	OSI-930 is a potent inhibitor of Kit, KDR and CSF-1R with IC50 of 80 nM, 9 nM and 15 nM, respectively...
Note	For research use only
Formula	$C_{22}H_{16}F_3N_3O_2S$
MW	443.44
CAS Number	[728033-96-3]
SMILES	<chem>O=C(C1=C(NCC2=CC=NC3=CC=CC=C23)C=CS1)NC4=CC=C(OC(F)(F)F)C=C4</chem>
Expiration Date	12 months from date of receipt.